

Exbivirumab ELISA Kit

Summary

Catalog No.	KDV03002
Alternative Names	19.79.5, HepeX-B (combination of exbivirumab and libivirumab), CAS: 569658-80-6
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	33.12 ng/mL
Range	78.125 - 5,000 ng/mL
Recovery	80-120%
Shipping	2-8 °C
Note	For Research Use Only.

Precision

CV

Data Image
