

Mabuterol/MAB ELISA Kit

Summary

Catalog No.	KSN56507
Alternative Names	MABUTEROL;Ambuterol;Benzenemethanol,4-amino-3-chloro- α -[[[(1,1-dimethylethyl)amino]methyl]-5-(trifluoromethyl)-(9CI);PB868Cl;4-Amino-3-chloro- α -[[[(1,1-dimethylethyl)amino]methyl]-5-(trifluoromethyl)benzChemicalbookenemethanol;1-[4-azanyl-3-chloro-5-(trifluoromethyl)phenyl]-2-(tert-butylamino)ethanol;4-Amino-3-chloro- α -[[[(tert-butyl)amino]methyl]-5-(trifluoromethyl)benzenemethanolhydrochloride
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Urine, Tissue, Feed
Assay type	Quantitative (Competitive)
Sensitivity	0.5 ng/mL
Range	1.56 - 100 ng/mL
Recovery	80-120%
Shipping	2-8 °C
Note	For research use only.

Precision

CV

Data Image
