

Anti-Chikungunya virus/CHIKV E2 hIgG ELISA Kit

Summary

Catalog No.	KAV09701
Alternative Names	Spike glycoprotein E2, E2, Anti-Chikungunya virus/CHIKV E2 hIgG
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Recovery	80-120%
Shipping	2-8 °C
Note	For research use only.

Precision

CV

Data Image
